Disposable Nitrile Examination Gloves

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Art. No.	Model	Size	Standard Package	Notified Body responsible for
23310320060	Worksafe Indigo Fine, powder free, blue	XS / 5-6	1/10 (150pcs/1500 pcs) 1/20 (50pcs/1500pcs)	EU Type Examination and Module B Centexbel 0493
23310320070	Worksafe Indigo Fine, powder free, blue	S / 6-7	1/10 (150pcs/1500 pcs) 1/20 (50pcs/1500pcs)	
23310320080	Worksafe Indigo Fine, powder free, blue	M / 7-8	1/10 (150pcs/1500 pcs) 1/20 (50pcs/1500pcs)	
23310320090	Worksafe Indigo Fine, powder free, blue	L / 8-9	1/10 (150pcs/1500 pcs) 1/20 (50pcs/1500pcs)	
23310320100	Worksafe Indigo Fine, powder free, blue	XL / 9-10	1/10 (150pcs/1500 pcs) 1/20 (50pcs/1500pcs)	

Notified Body which carried out:
EU Type Examination Certificate According to Module B Issued by:
Centexbel, Technologiepark 70, 9052, Zwijnaarde, Belgium (Notified Body No. 0493)
The Module D is under Supervision of:
SGS Fimko Oy P.O Box 30 (Sarkiniementie 3) 00211, Helsinki, Finland (Notified body No. 0598)

The product is manufactured by
Procurator AB, P.O. Box 9504, SE-200 39 Malmö, Sweden
Tel +46-(0) 10 60 40 000, www.procurator.com

EU Declaration of Conformity can be downloaded from <http://doc.worksafe.com>

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FR

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Made in Malaysia

EN ISO 374-1:2016 - Permeation levels are based on breakthrough times


EN 16523-1:2015 - Resistance to chemical permeation

EN ISO 374-4:2013 - Resistance to chemical degradation

EN ISO 374-5:2016 - Protection against bacteria and fungi, viruses

* Fig. 1		Breakthrough Performance level	Achieved breakthrough time (min)
ISO 374-1:2016  ABCDEF GHIJ KLMNOPST	A. Methanol	J. n-Heptane	>10min
	B. Acetone	K. Sodium hydroxide 40%	>30min
	C. Acetonitrile	L. Sulphuric acid 96%	>60min
	D. Dichloromethane	M. Nitric acid 65%	>120min
	E. Carbon disulphide	N. Acetic acid 99%	>240min
	F. Toluene	O. Ammonium hydroxide 25%	>480min
G. Diethylamine		P. Hydrogen peroxide 30%	
H. Tetrahydrofurane		S. Hydrofluoric acid 40%	
I. Ethyl acetate		T. Formaldehyde 37%	

REF/ART NO.	Resistance to Permeation by chemicals	Breakthrough Performance Level (BS EN 16523-1:2015)	Mean Degradation(%) (BS EN 374-4:2013)	ISO 374-5:2016  VIRUS
2331032 Worksafe Indigo Fine, powder free, blue	(U) n-Heptane	2	50,4	Protection against bacteria and fungi - Pass
	K. Sodium Hydroxide 40%	6	6,7	Protection against viruses - Pass
	(T) 37% Formaldehyde	6	38,5	Performance Level Class 2 (AQL 1.5, G1)

EN ISO 21420:2020, glove sizing:

Choosing the right size of gloves is key as it affects the dexterity and comfort of the working hand. Please use the above size chart.

GB
The glove is in compliance with MDR (EU) 2017/745 as a Class I medical device. The gloves are tested in compliance with EN455 -1, 2, 3 & 4, Disposable gloves for medical use. EN455-1 Medical gloves for single use - Part 1:Requirement and testing for freedom from holes. EN455-2 Medical gloves for single use - Part 2: Requirements and testing for physical properties. EN455-3 Medical glove for single use - Part 3: Requirements and testing for biological evaluation. EN455-4 Medical glove for single use - Part 4: Requirements and testing for shelf life determination.

AQL= Acceptable quality level that indicates the percentage of leaky gloves per production run (or quantity).

PERFORMANCE LEVEL	AQL	INSPECTION LEVEL
Level 3	<0,65	G1
Level 2	<1,5	G1
Level 1	<4,0	S4

EN ISO 374-1:2016
Protection against chemicals and microorganisms.

ISO 374-1:2016
Type A/B/C

ABCDEF GHIJ KLMNOPST

EN ISO 374-5:2016
Protection against bacteria and fungi. Protection against viruses.

ISO 374-5:2016

VIRUS

**Fig 1*

ISO 374-1:2016
Type C

ABCDEF GHIJ KLMNOPST

ISO 374-1:2016
Type B

ABCDEF GHIJ KLMNOPST

ISO 374-1:2016
Type A

ABCDEF GHIJ KLMNOPST

Type C - At least Breakthrough Performance Level Class 1 (more than 10 minutes) against at least 1 chemical on the list. **Fig 1*

Type B - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 3 chemicals on the list. **Fig 1*

Type A - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 6 chemicals on the list. **Fig 1*

The gloves comply with PPE Regulation 2016/425 on personal protective equipment, Category III. The gloves are tested in accordance with EN ISO 21420:2020, EN420:2003+A1:2009 General requirements for gloves, EN ISO 374-1:2016, EN16523-1:2015 and EN ISO 374-5:2016.
The glove is intended for handling foodstuffs and cleaning, and to protect hands from dirt, minimal risks. "Fit for special purpose gloves" because they are to be used to protect the hand only from chemical splashes when handling chemicals. Do not use these gloves when protection in the cuff area is needed. Due to this, the glove does not fully correspond with the standard lengths declared in EN 21420. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and

degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
EN 374-4 Degradation: Levels indicate the change in puncture resistance of the glove after exposure to the tested chemical. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Water vapor transmission: if the glove is not breathable, and therefore not tested, you can improve the comfort by limiting time of use. Dexterity: Level 5. Donning: Ensure the selection of appropriate glove size, check for any damages or contamination prior to each use. Discard if there is any damage or contaminant. Doffing: remove gloves without the contaminant contacting bare hand or skin. Before usage, inspect the gloves for any defects or imperfections. Protect from heat, moist, light and ozone. Keep unused gloves in the original packaging. To be stored in a cool dark room. Disposable glove must be discarded after use. The product material is not known to cause allergic reactions. If irritation of the skin should occur, wash the affected area with mild soap and water. Seek medical attention if the irritation persists.

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