WORKSAFE

Disposable Nitrile Examination Gloves

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AQL 1,5 EN 455 - 1, 2, 3, 4

These gloves are tested and approved by Notified Body: 2777 module B
SATRA Technology Europe Ltd Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland

Notified Body responsible for Module D Conformity Assessment:
Satra Technology Europe Ltd. Notified Body No. 2777 Bracetown Business Park, Clonee, Dublin 15 Dublin, Ireland

EN ISO 374-1:2016+A1:2018 fig 1

A. Methanol 67-56-1
B. Acetone 67-64-1
C. Acetonitrile 75-05-8
D. Dichloromethane 75-09-2
E. Carbon disulphide 75-15-0
F. Toluene 108-88-3
G. Diethylamine 109-89-7
H. Tetrahydrofurane 109-99-9
I. Ethyl acetate 141-78-6
J. n-Heptane 142-82-5
K. Sodium hydroxide 40% 1310-73-2
L. Sulphuric acid 96% 7664-93-9
M. Nitric acid 65% 7697-37-2
N. Acetic acid 99% 64-19-7
O. Ammonium hydroxide 25% 1336-21-6
P. Hydrogen peroxide 30% 7722-84-1
S. Hydrofluoric acid 40% 7664-39-3
T. Formaldehyde 37% 50-00-0

Breakthrough Performance level	Achieved breakthrough time (min)
Class 1	>10min
Class 2	>30min
Class 3	>60min
Class 4	>120min
Class 5	>240min
Class 6	>480min

AQL Fig. 2

Performance level	AQL	Inspection level
Level 3	<0,65	G1
Level 2	<1,5	G1
Level 1	<4,0	S4

Art. No.	EN ISO 374-1:2016+A1:2018 Permeation levels are based on breakthrough times	Breakthrough Performance Level	EN ISO 374-4:2019 Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical Resistance to chemical degradation (%)	EN ISO 374-5:2016 VIRUS
2155897, 2155898, 2155899, 2155900, 2155901, 2162247, 2162248, 2162249, 2162250, 2162251,	40% Sodium Hydroxide (K) 30% Hydrogen Peroxide (P) 37% Formaldehyde (T)	6 6 6	-95,4 26,5 6,1	Protection against bacteria and fungi - Pass Protection against viruses - Pass

GB

EN ISO 374-1:2016+A1:2018

Type A

UVWXYZ

Type A - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 6 chemicals on the list. *Fig 1

EN ISO 374-1:2016+A1:2018

Type B

XYZ

Type B - At least Breakthrough Performance Level Class 2 (more than 30 minutes) against at least 3 chemicals on the list. *Fig 1

EN ISO 374-1:2016+A1:2018

Type C

X

Type C - At least Breakthrough Performance Level Class 1 (more than 10 minutes) against at least 1 chemical on the list. *Fig 1

EN ISO 374-1:2016+A1:2018 Protection against chemicals and microorganisms.

EN ISO 374-5:2016 Protection against bacteria and fungi. Protection against viruses.

AQL= Acceptable quality level that indicates the percentage of leaky gloves per production run (or quantity).

This product complies with the requirements of the Regulation (EU) 2016/425 on Personal Protective Equipment Category III, MDR (EU) 2017/745 concerning Medical Devices Class I, Regulation (EC) No 1935/2004 on Materials and Articles intended to come into Contact with Food.

The gloves are tested in accordance with EN ISO 21420:2020 General requirements for gloves, EN ISO 374-1:2016+A1:2018 and EN ISO 374-5:2016. EN 455-1, EN 455-2, EN 455-3, EN 455-4

Fit for special purpose gloves" because they are to be used to protect the hand only from chemical splashes when handling chemicals. Do not use these gloves

when protection in the cuff area is needed. Due to this, the glove does not fully correspond with the standard lengths declared in EN ISO21420:2020. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Before usage, inspect the gloves for any defects or imperfections. Protect from heat, moist, light and ozone. Keep unused gloves in the original packaging. To be stored in a cool dark room. Disposable glove must be discarded after use. The product material is not known to cause allergic reactions. If irritation of the skin should occur, wash the affected area with mild soap and water. Seek medical attention if the irritation persists.

Donning: 1. Remove all hand and wrist jewelry, and wash the hands before donning. 2. The user only touch the inside of the gloves when putting on the first glove. 4. The wearer uses the gloved hand to slip the other glove onto the other hand. 5. Once both gloves are on, the users can touch the outside of the gloves to ensure a proper fit.

Doffing: 1. Start by grabbing the outside of the glove on one hand on the palm side near the cuff. 2. Pull the glove off the hand and place it in the gloved hand, balling it up. 3. Slip two fingers under the cuff of the other hand glove and carefully peel it off the hand turning the remaining glove inside out as it is removed and in turn encasing the first glove. 4. The gloves can be disposed.

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